

# **dPette+** Multi-functional 8-channel Electronic Pipette

DLAB multi-functional high performance 8-channel electronic pipette offers productive pipetting with its easy to understand operations. Its uniquely lightweight streamlined design ensures an effortless transfer of multiple samples with an increased throughput and data reproducibility.



#### Three easy steps to operate dPette+

1

Long press the Parameter Knob for 2 seconds to Start

Rotate the knob back and forth quickly to Switch Pipetting,Continuous Dispensing, and other Function Settings

2

3

Quickly turn the Parameter Knob to unlock it, turn it to adjust the parameter such as volume, press it to pipette, long press it to mix

### j iquid Handling

#### **Pipetting function**

Parameter Knob press → Pipetting, Long press → Mixing

# PI ST ST Pipetting volume I I I I I Aspirate I Dispense

#### **Stepper function**

PI

Maximum steps = Pipette nominal volume / One step dispensing volume

**One-step** 

dispensing volume

nle numbers

100 µL×3

#### **Other function**

Pipetting speed adjustment Key tone ON/OFF



**Contact Charging** 

# USB Charging

**Supports dual charging mode** 

#### **Support automatic Calibration**

Connect it to the computer, with the free software provided for calibration

## **Specifications**



Complimentary Pipette Holder

| Channels | Volume Range | <b>increment</b><br>μL | <b>Test Volume</b><br>μL | Systematic Error |       | Random Error |       |
|----------|--------------|------------------------|--------------------------|------------------|-------|--------------|-------|
|          | μL           |                        |                          | μL               | %     | μL           | %     |
| 8        | 0.5-10       | 0.01                   | 10                       | ±0.20            | ±2.00 | ±0.10        | ±1.00 |
|          |              |                        | 5                        | ±0.20            | ±4.00 | ±0.10        | ±2.00 |
|          |              |                        | 1                        | ±0.08            | ±8.00 | ±0.05        | ±5.00 |
| 8        | 10-100       | 0.1                    | 100                      | ±0.80            | ±0.80 | ±0.30        | ±0.30 |
|          |              |                        | 50                       | ±0.50            | ±1.00 | ±0.40        | ±0.80 |
|          |              |                        | 10                       | ±0.30            | ±3.00 | ±0.20        | ±2.00 |
| 8        | 15-300       | 1                      | 300                      | ±1.80            | ±0.60 | ±0.90        | ±0.30 |
|          |              |                        | 150                      | ±1.50            | ±1.00 | ±0.75        | ±0.50 |
|          |              |                        | 15                       | ±0.60            | ±4.00 | ±0.15        | ±1.00 |
| 8        | 30-300       | 1                      | 300                      | ±1.80            | ±0.60 | ±0.90        | ±0.30 |
|          |              |                        | 150                      | ±1.50            | ±1.00 | ±0.75        | ±0.50 |
|          |              |                        | 30                       | ±0.90            | ±3.00 | ±0.30        | ±1.00 |

\* DLAB specifications are used as guidelines and the user calibration should refer to the industrial standard ISO 8655.