

# Optional Add-ons

## Distillation System Connection (For RE100-Pro&RE100-S)



C410



### C410 Chemistry-design diaphragm pump








#### Specifications

Maximum Vacuum:13mbar  
 Work Mode: Continuous Operation  
 Nozzle Size:8mm, Power: 95W  
 Motor speed:1450rpm  
 Flow Rate:

mbar	1000	880	750	600	480	280
L/min	25	22	20	12	7.5	5

Noise:50dB

### Seals and other accessories

Cat. No.	Description	
18204186	vacuum seal ( PTFE+Viton)	
18204187	Evaporation bottle spring Used with 18204186	
18900190	Evaporation flask clip(red) NS29/42	
18900188	Evaporation flask clip(green) NS24/40	
18901001	Evaporation flask clip(Black) , NS24/40, Adapted to 3L/5L Evaporating flask	
18901002	Evaporation flask clip(Black) , NS29/32, Adapted to 3L/5L Evaporating flask	
18900871	Cooling water silicone hose line	
18900879	Secondary Condenser Component	
18900744	Rotary Evaporator Protective Shield Assembly (compatible with RE100-Pro)	
18900789	Rotary Evaporator Protective Shield (RE100-S)	

### Regulator /Moisture

Cat. No.	Description	
18900186	Material: PVDF, PTFE, ABS, Borosilicate Glass, Stainless Steel	