

SK-R330-Pro

**Large Multifunctional
Transference Rocking Shaker**

SK-R1807-S

**Intelligent Transference
Decolorizing Shaker**

Features

- Rocking motion angle 9°(R330-Pro) and 7°(R1807-S)
- Two options:
10kg (R330-Pro) and speed range 10-70rpm,
2kg (R1807-S) and speed range 10-80rpm
- Double LCDs/LEDs for independent display of time and speed
- Continuous or timer operation mode
- Electronic time switching clock: 1min-99h59min
- Wide range of platforms for use with a variety of vessels(R330-Pro)
- Fixed work plate(R1807-S)
- Simple appearance, save operating space

Applications

Designed for time-dependent and low rpm sensitive shaking applications such as ELISA washes, enzyme immunoassays, protein synthesis, hybridization staining/destaining gels, immune precipitation and Western blots, cell culture, pharmaceutical profiling, etc.





Rocking Motion
















SK-R330-Pro











SK-R1807-S

**SK-R1807-S (fixed plate)****SK-R1807-S (multiple plate are optional)**

Specifications	SK-R330-Pro	SK-R1807-S
Tilt angle	9°	7°
Max. load capacity(with platform)	10kg	2kg
Motor type	DC motor	DC motor
Motor input power	40W	20W
Motor output power	24W	16W
Speed range	10-70rpm	10-80rpm
Speed display	LCD	LED
Timer display	LCD	LED
Time settings range	1min-99h59min	1min-99h59min
Voltage,Frequency	100-240V, 50/60Hz	100-240V, 50/60Hz
Power	50W	20W
Weight	9kg	3.1kg
Permissible ambient temperature and humidity	5-40°C,80%RH	5-40°C,80%RH
Dimension[W×D×H]	360×410×200mm	294×348×123mm
Protection class	IP21	IP21

Cat. No.	Description	SK-0330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-0330-M	SK-0180-Pro	SK-L180-Pro	SK-0180-S (Multifunc)	SK-R1807-S (Multifunc)
 18900155	SK330.5, Dish platform with non-slip mat and fixed string (8pcs) Dimension[WxD]: 30x28.5cm	■	■	■	■					
 18202835	Silicone mattress for 16× 5-15mL tubes, use with SK330.5	■	■	■	■					
 18202846	Silicone mattress for 6 × 50mL tubes, use with SK330.5	■	■	■	■					
 18900185	SK180.4, Dish attachment with non-slip mat Dimension[WxD]: 24x21.5cm						■	■	■	■
 18900869	SK330.4 Dish platform with non-slip mat and fixed string(8 pcs), 34X34cm					■				
 18900859	SK330.6 Transparent shield, serologic test for syphilis, 41X39.5X10cm	■	■			■				
 18900860	SK180.5 Transparent Cover for Syphilis Serum Test, 32X32X8.3cm						■	■		
 18202852	Silicone mattress for 10× 5-10mL tubes						■	■	■	■
 18202853	Silicone mattress for 4 × 50mL tubes						■	■	■	■
 18900027	SK330.1, Universal attachment with 4 bars Size: 32x32cm Capacity for flask: 25mLx21pcs, 50mLx15pcs, 100mLx12pcs, 250mLx6pcs, 500x4pcs	■	■	■	■					
 18900040	SK330.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*3pcs 250mL*4pcs 100/125mL*4pcs 50mL*5pcs	■	■	■	■					
 18900025	SK180.1, Universal attachment with 3 bars Size: 25x24.5cm Capacity for flask: 25mLx10pcs, 50mLx8pcs, 100mLx6pcs, 250mLx4pcs, 500x2pcs						■	■	■	■
 18900038	SK180.3 Lengthways roller attachment, Mainly used for separation funnel oscillation extraction applications 500mL*2pcs 250mL*3pcs 100/125mL*3pcs 50mL*4pcs						■	■	■	■

Cat. No.	Description	SK-O330-Pro	SK-L330-Pro	SK-R330-Pro	SK-D3309-Pro	SK-O330-M	SK-O180-Pro	SK-L180-Pro	SK-O180-S (Multifunc)	SK-R1807-S (Multifunc)
 18900036	Bar for SK330.1 universal attachment	■	■	■	■					
 18900037	Bar for SK180.1 universal attachment						■	■	■	■
 18900280	Connect sticks for double dishes shaker	■	■	■	■	■	■	■	■	■
 18900028	SK330.2, Fixing clip attachment Size: 32 x 32cm	■	■	■	■					
 18900026	SK180.2, Fixing clip attachment Size: 25 x 24cm						■	■	■	■
 18900029	Fixing clip for flask volume 25mL, used with SK330.2[28pcs] SK180.2[9pcs]	■	■	■	■		■	■	■	■
 18900030	Fixing clip for flask volume 50mL, used with SK330.2[16pcs] SK180.2[9pcs]	■	■	■	■		■	■	■	■
 18900031	Fixing clip for flask volume 100mL, used with SK330.2[16pcs] SK180.2 [9pcs]	■	■	■	■		■	■	■	■
 18900032	Fixing clip for flask volume 200/250mL, used with SK330.2[9pcs] SK180.2[4pcs]	■	■	■	■		■	■	■	■
 18900033	Fixing clip for flask volume 500mL, used with SK330.2[9pcs] SK180.2[2pcs]	■	■	■	■		■	■	■	■