

# dTrite

## Electronic Burette

Digital Burette - dTrite provides accurate, precise and convenient bottle-top titration with optimum operator safety. It combines both magnetic stirrer and Burette functions, making it easier for the user to achieve reproducible results in the field of analytical chemistry, food industry, water analysis etc.

## Features

- High quality motor with excellent precision and accuracy
- Electronic control provides fatigue-free operation
- Remote control panel prevents manual disturbance during operation
- Built-in magnetic stirrer for a complete and easy operation
- Large LCD display parameters



## Specifications

Volume Range	0.01mL-99.99mL, Max piston lift is 10mL, resolution is 10 $\mu$ L
Volume Accuracy	R=0.2% CV=0.07%
Speed	16 Levels
Operating Temperature Range	10 $^{\circ}$ C-30 $^{\circ}$ C
Quality Control	DIN EN ISO 8655
Control Type	External control, including the stirrer and the burette program

\* Notice: R: Relative error; CV: coefficient of variation.