

Mini Pipette

Mini Pipette is specially customized for fast pipetting with fixed volume and it is the first choice for IVD assay kit accessory. Ultra-light body, precise pipetting, smooth aspiration and dispensing bring an easier pipetting experience to customers.

- Ergonomic design, compact body, with easy operation;
- Easy to press, consistent with standard manual pipettes in operating habits;
- Multiple optional volume from 10 to 300ul, marked with different colors;
- Compatible with the standard 200ul/300ul tips;
- Accurate pipetting and reliable quality;
- Support customization;



Standard style



New style

Color	Light Green	Grey	Dark Green	Yellow	Orange	Gray	Red	Blue	Grayish green	Other
Mini Pipette	10μL	20μL	25μL	30μL	50μL	60μL	75μL	100μL	200μL	OEM
Color	Dark Green	Orange	Red	Light Blue	Blue	Gray	Purple	Other	—	—
New Mini Pipette	25	50	75	80	100	200	300	OEM	—	—

Levo E Pipette Pump

Designed for single-handed operation, Pipette Pump identifies itself with a thumbwheel, allowing precise aspiration or dispensing. The feature of side lever enables complete liquid dispensing, without a droplet remaining.



Features

- Volume capacity 2mL, 5mL, 10mL and 25mL
- Colour-coded by volume with Blue, yellow, green and red
- Thumb wheel guarantees precise operation
- Resistant to acids and alkalis
- Easy for cleaning and maintenance