

# MS-T-S15

## 15-Position magnetic stirrer

### Features

- LCD display for precise speed and time control.
- Maintenance free DC brushless motor for long life and high stirring power.
- Aluminum laser engraving work plate surface, anti-abrasive, skid and corrosion proof
- Strong magnets ensure fast and intense stirring
- 15 stirring positions are divided into 3 independent groups and each group can be controlled respectively at different speeds and with different timer settings.
- Buzzer function : When the time countdown ends, the buzzer will remind the user by sound.



Specifications	MS-T-S15
Work plate Dimension [WXD]	460x304mm
Work plate material	Aluminum
Motor type	DC brushless
Stirring positions	15
Max. stirring quantity [HzO]	0.4Lx15
Max. magnetic bar[length]	40mm
Speed control accuracy	±10rpm
Speed、Time display	LCD
Speed range	100-1500rpm
Timer function	1min-99h59min
Protection class	IP21
Power	54W
Voltage、Frequency	100-240V, 50Hz/60Hz
Dimension [DxWxH] [mm]	615x307x90mm
Weight	9.2kg
Permissible ambient temperature and humidity	5-40°C, 80%RH



0.4Lx15



100-1500rpm



Timer



Skid resistance



Multiple position



Brushless DC motor